

Empowering your research with services from the CCR Genomics Core

Madeline Wong, Ph.D.
Staff Scientist (Office of Science Technology and Research, OSTR)
Building 41, Room D310

Disclaimer


We did not receive any financial support from any suppliers or manufacturers mentioned and do not have any conflict of interest

Background

- 3 NCI CCR Genomics/ sequencing cores
 - Sequencing Facility
(Bao Tran, Frederick)
 - Genomics Technology
Laboratory
(Xiaolin Wu, Frederick)
 - Genomics Core
(Liz Conner, Bethesda)

Background

- "Open access" facility with goal of providing efficiency and quality of a centralized facility with speed and convenience
- Help identify, evaluate and make available new technologies as they emerge
- Provide services to investigators in NCI CCR scientific community, the Clinical Center and other NIH institutes including NIAID, NHGRI, NHLBI, NINDS, NEI, NIA, NIMH, NIDDK, and NICHD



Staff & Location

Bldg. 41, Room D310



Collaborative Protein
Technology Resource
(CPTR)

Single Cell Analysis Facility
(SCAF)


Genomics Core
(GC)

Core members:

- Liz Conner, Manager
- Steve Shema
- Qin Wei
- Madeline Wong
- Desiree Tillo,
Bioinformatics
Consultant

Facility Usage and Access

Investigators can register and submit samples through our website at: <https://nci.corefacilities.org/account/login>



- Users complete iLAB request for Sanger samples
- Other applications the Core Staff will enter the iLAB request
- Billing every quarter of each fiscal year (DEC, APR & AUG) via CAN transfer
- NO POTS order and any other type of approval

Services

- **Analytical & Preparative Electrophoresis**
- Droplet Digital PCR (ddPCR)
- Digital Gene Expression (Nanostring nCounter)
- Sanger Sequencing
- Next-Generation Sequencing (NGS) – illumina platforms
- Digital Spatial Profiling (NanoString GeoMx DSP)
- Next-Generation Sequencing (NGS) – Oxford Nanopore

Service: Analytical & Preparative Electrophoresis

 sage science



PippinHT

High Throughput Sizing

- Up to 24 samples in one run
- 3% agarose, 100-250 bp
- 2% agarose, 100-600 bp
- 1.5% agarose, 300-1.5 kb
- 0.75% agarose, High-pass size selection

• **Cost recovery service (\$10.00/sample) with same day sample delivery**

 Agilent Technologies



4150 & 4200 Tape Station

Automated Electrophoresis for DNA & RNA samples to analyze size, quantity & integrity (RNA RIN number)

- Genomic DNA Tape
- D1000 DNA standard and high sensitivity
- D5000 DNA standard and high sensitivity
- RNA standard and high sensitivity

• **Cost recovery service (\$5 - 6.00/sample) with same day data delivery**

Services

- Analytical & Preparative Electrophoresis
- **Droplet Digital PCR (ddPCR)**
- Digital Gene Expression (Nanostring nCounter)
- Sanger Sequencing
- Next-Generation Sequencing (NGS) – illumina platforms
- Digital Spatial Profiling (NanoString GeoMx DSP)
- Next-Generation Sequencing (NGS) – Oxford Nanopore

Service : Droplet Digital PCR (ddPCR)

BIO-RAD QX200 Droplet Digital PCR System




- Provides absolute quantification of target molecules without the use of standard curves
- Chemistries
 - EvaGreen
 - Probe
- Common Applications
 - Expression
 - Mutation
 - Copy Number Variation

Service : Droplet Digital PCR (ddPCR)

BIO-RAD QX200 Droplet Digital PCR System

- Users provide DNA template with primers and Probe (FAM & HEX Probe Assay: 2 genes) or DNA template with primers (EvaGreen Assay)
- Core Staff adds the Supermix, Generate droplets, Perform PCR and Read
- Biorad Analysis software provided free to Users
- Users will analyze their results with Biorad support if needed



• Cost recovery service (\$7.00/sample in batch of 8 samples) with same day or next day data delivery

Service : Droplet Digital PCR (ddPCR)

BIO-RAD QX200 Droplet Digital PCR System

- Probe Assay – 2 colors systems (FAM & HEX)

▶ **COMING SOON**



BIO-RAD QX600 : 6 plex (colors) system

Increase power of experiment – multiplex to detect > 2 genes of interest in a single well

Services

- Analytical & Preparative Electrophoresis
- Droplet Digital PCR (ddPCR)
- **Digital Gene Expression (Nanostring nCounter)**
- Sanger Sequencing
- Next-Generation Sequencing (NGS) – illumina platforms
- Digital Spatial Profiling (NanoString GeoMx DSP)
- Next-Generation Sequencing (NGS) – Oxford Nanopore

Service : NanoString Digital Gene Expression

nanoString iCounter Analysis System



- Novel digital color-coded barcode technology (direct multiplexed measurement of gene expression)
- High levels of precision and sensitivity
- Users contact NanoString directly to purchase or design a CodeSet tailored to their experiment (Pre-made panels are available for Gene Expression, Gene Fusions, miRNA, CNV, SNV and Protein)
- Samples are processed after hybridization in the Core



1 Hybridize

- Flexible sample requirements
- Only 4 pipetting steps
- No Amplification
- 800 hybridizations in a single tube



2 Purify



3 Count


- Sensitive
- Precise
- Quantitative
- Simple

• **Cost recovery service (\$24.00 / sample in batch of 12 samples) with results within 2 business days from the start of hybridization**

Services

- Analytical & Preparative Electrophoresis
- Droplet Digital PCR (ddPCR)
- Digital Gene Expression (Nanostring nCounter)
- **Sanger Sequencing**
- Next-Generation Sequencing (NGS) – illumina platforms
- Digital Spatial Profiling (NanoString GeoMx DSP)
- Next-Generation Sequencing (NGS) – Oxford Nanopore

Service : Sanger Sequencing



- Applied Biosystems DNA Analyzer **ThermoFisher**
▪ 24-capillary 3500XL (2 units) **SCIENTIFIC**
- Standard plasmid and PCR products as well as bacterial genomic DNA, BACs and cosmids. Fragment analysis – analyze amplified PCR products
- Investigators provide us with DNA template + primer and we perform the reaction and electrophoresis or customers do their own reactions for electrophoresis only
- Data is analyzed and placed in a secure folder on the NIH network and investigators are notified via email when their data becomes available

• **Cost recovery service (\$2 - 6.00/sample) & typically 1 business day data delivery for samples received by 1 pm**

Services

- Analytical & Preparative Electrophoresis
- Droplet Digital PCR (ddPCR)
- Digital Gene Expression (Nanostring nCounter)
- Sanger Sequencing
- **Next-Generation Sequencing (NGS) – illumina platforms**
- Digital Spatial Profiling (NanoString GeoMx DSP)
- Next-Generation Sequencing (NGS) – Oxford Nanopore

Service : Next-Generation Sequencing (NGS)

illumina platforms



iSeq

- 1 kit
- 5 M max clusters
- Up to 2 x 150 bp



MiSeq

- 8 different kits available
- 25 M max cluster
- Up to 2 x 300 bp



NextSeq550

- 5 different kits available
- 400 M max clusters
- Up to 2 x 150 bp



NextSeq2000

- 11 different kits available
- 1.1 B max cluster
- Up to 2 x 300 bp

• **Cost recovery service (\$400 – 5700/kit) & data delivery < 1 week**

- Library provided by Users (RNASeq, CHIPSeq, NanoString DSP, 10X Genomics – 3', 5', Immune and Visium, Amplicon or Custom Library prep)
- Library prepared by us (RNASeq, CHIPSeq, Bacterial genome, NanoString DSP)
- Different Users' samples are not

Service : Next-Generation Sequencing (NGS)

illumina platforms

- For >1.1B reads, Users are advised to discuss with the Sequencing Facility in Frederick for production scale NovaSeq run



NovaSeq6000

20B max cluster

Up to 2 x250 bp



NovaSeq X

52B max cluster


Up to 2 x150 bp

Services

- Analytical & Preparative Electrophoresis
- Droplet Digital PCR (ddPCR)
- Digital Gene Expression (Nanostring nCounter)
- Sanger Sequencing
- Next-Generation Sequencing (NGS) – illumina platforms
- **Digital Spatial Profiling (NanoString GeoMx DSP)**
- Next-Generation Sequencing (NGS) – Oxford Nanopore

Service : GeoMx Digital Spatial Profiler

- Shared between Collaborative Protein Technology Resource (CPTer) and Genomics Core based on expertise
- High-plex spatial profiling to assess protein and/or RNA within tissue samples from a single FFPE or



- Interested Users would contact CPTer or Genomics Core to obtain more information

Assay types:

(DSP)

PROTEIN

- Protein with nCounter readout
 - Human Immuno oncology
 - Human Neuroscience
 - Mouse Immuno oncology
 - Mouse Neuroscience
- Protein with NGS readout (illumina platforms)
 - Human core + immuno, oncology or neuro panels

RNA

- RNA with nCounter readout: 94 transcripts (73 core panel genes + 5 housekeeping + 6 neg control + 10 custom)
 - Human Immune pathways
- RNA with NGS readout (illumina platforms)
 - Human Cancer Transcriptome Atlas (hu CTA): ~ 1800 transcripts (~5 probes per gene)
 - Human Whole Transcriptome Atlas (hu WTA): ~ 18,000 transcripts (one probe per gene)
 - Mouse Whole Transcriptome Atlas (m WTA): ~18,000 transcripts (one probe per gene)

Services

- Analytical & Preparative Electrophoresis
- Droplet Digital PCR (ddPCR)
- Digital Gene Expression (Nanostring nCounter)
- Sanger Sequencing
- Next-Generation Sequencing (NGS) – illumina platforms
- Digital Spatial Profiling (NanoString GeoMx DSP)
- **Next-Generation Sequencing (NGS) – Oxford Nanopore**

Service : Next Generation Sequencing (NGS)

Long-read sequencing



PacBio



Sequencing Facility
in Frederick

SEQUEL IIe SYSTEM

Oxford
NANOPORE
Technologies



Sequencing Facility in
Frederick
Genomics Core

Service : Next Generation Sequencing (NGS)

Oxford **NANOPORE** Technologies Long-read sequencing



	Flongle	MinION	GridION	P2	PromethION
WGS - small genomes	Yes: Low to medium plex	Yes: Low to medium plex	Yes: Low to medium plex	Yes: Highly multiplexed	Yes: Highly multiplexed
WGS - large genomes	Low pass	Low pass	Yes	Yes	Yes
Targeted sequencing	Yes: Low to medium plex, and adaptive sampling			Yes: Highly multiplexed and adaptive sampling	
Metagenomics	Yes: Quantitative species ID			Yes: Quantitative species ID	
RNA sequencing	Yes: Isoform & expression from same experiment			Yes: Isoform & expression from same experiment, including low abundance transcripts and single cell	
Epigenetics	Yes - when native DNA used			Yes - when native DNA used	

Service : Next Generation Sequencing (NGS)



Long-read sequencing

Kits

FLO-MIN106D flow cells are suitable for most sequencing kits:

- Ligation Sequencing Kit (SQK-LSK110)
- Ligation Sequencing Kit (SQK-LSK109)
- Cas9 Sequencing Kit (SQK-CS9109)
- PCR-cDNA Sequencing Kit (SQK-PCS111)
- PCR-cDNA Sequencing Kit (SQK-PCS109)
- PCR-cDNA Barcoding Kit (SQK-PCB109)
- Direct cDNA Sequencing Kit (SQK-DCS109)
- Direct RNA Sequencing Kit (SQK-RNA002)
- Rapid Sequencing Kit (SQK-RAD004)
- Rapid Barcoding Kit (SQK-RBK004)
- Rapid PCR Barcoding Kit (SQK-RPB004)
- 16S Barcoding Kit (SQK-RAB204)
- PCR Sequencing Kit (SQK-PSK004)
- PCR Barcoding Kit (SQK-PBK004)
- Field Sequencing Kit (SQK-LRK001)
- Ligation Sequencing Kit (SQK-LSK112) - *legacy*
- Ligation Sequencing Kit XL (SQK-LSK112-XL) - *legacy*
- Native Barcoding Kit 24 (SQK-NBD112.24) - *legacy*
- Native Barcoding Kit 96 (SQK-NBD112.96) - *legacy*

N⁶-methyladenosine (m⁶A) mRNA modification

A specific gene expression

Whole plasmid sequencing

Viral genome sequencing

Software

Basecalling:

- MinKNOW
- Guppy

User analyzing the data

as .fast5 and FASTQ files.

Downstream analysis:

• EPI2ME

The gene was not detected

tools and pipelines

customer-developed tools and pipelines

User analyzing the data

FLO-MIN114 flow cells can only be used with the V14 Sequencing Kit:


- Ligation Sequencing Kit V14 (SQK-LSK114)

Service : Next Generation Sequencing (NGS)



Long-read sequencing

- Rapid PCR Barcoding Kit (SQK-RPB004) ← **Whole plasmid sequencing**



- Advantage to the CCR community if we offer the service – no POTS order or any other approval – will be similar to requesting Sanger sequencing
- Maybe useful for difficult to assemble plasmids



Affordable, fast, **Whole Plasmid Sequencing** using Oxford Nanopore Technologies is now available for **\$15 per sample**. We offer **free** sample pick-up in DC area or you may ship directly to our facility, and **NO** sequencing primers are required!

Service : Next Generation Sequencing (NGS)





• Rapid PCR Barcoding Kit (SQK-RPB004)

← **Whole plasmid sequencing**

Sample: p304_ABE_E59A_SpG

• obtained single assembly

Completed successfully



Service : Next Generation Sequencing (NGS)



Long-read sequencing

- Rapid PCR Barcoding Kit (SQK-RPB004) ← Whole plasmid sequencing, amplicon or bacterial

plasmidsaurus
Affordable, fast, long-read sequencing
Certified by Oxford Nanopore Technologies

Whole Plasmid Sequencing for **\$15**

- Complete, accurate plasmid sequencing from 2.5 kb to 300 kb
- Each plasmid is fully sequenced by long reads to high consensus accuracy
- Results returned within one day

Learn more
Order

plasmidsaurus
Affordable, fast, long-read sequencing
Certified by Oxford Nanopore Technologies

Bacterial Genome	Linear/Amplicon
We sequence bacterial genomes up to 12 Mb with long sequencing reads and return an assembly with annotations. From \$90/sample , ~5 days	Long-read linear/amplicon sequencing from 600 bp to 125 kb . From \$15/sample , one day

Affordable, fast, **Whole Plasmid Sequencing** using Oxford Nanopore Technologies is now available for **\$15 per sample**. We offer **free** sample pick-up in DC area or you may ship directly to our facility, and **NO** sequencing primers are required!


- The cost – we are unable to match the price unless we have enough samples to

Service : Next Generation Sequencing (NGS)



Long-read sequencing

- Rapid PCR Barcoding Kit (SQK-RPB004) ← **Whole plasmid sequencing**
 - Currently comparing the Sanger, short-read illumina and long-read Oxford Nanopore data
 - higher error rate, especially around homopolymers (illumina generally comparable to Sanger)
 - Need more testing



140bp gap in ONT assembly

Service : Next Generation Sequencing (NGS)



Long-read sequencing



- Hosted a seminar on Feb 15, 2023: Single cell sequencing with Oxford Nanopore and 10X Genomics
- collaboration with SCAF
- Other applications

	Flongle	MinION	GridION	P2	PromethION
WGS - small genomes		Yes: Low to medium plex	Yes: Low to medium plex	Yes: Low to medium plex	Yes: Highly multiplexed
WGS - large genomes		Low pass	Low pass	Yes	Yes
Targeted sequencing		Yes: Low to medium plex, and adaptive sampling			Yes: Highly multiplexed and adaptive sampling
Metagenomics		Yes: Quantitative species ID			Yes: Quantitative species ID
RNA sequencing		Yes: Isoform & expression from same experiment			Yes: Isoform & expression from same experiment, including low abundance transcripts and single cell
Epigenetics		Yes - when native DNA used			Yes - when native DNA used

Service : Next Generation Sequencing (NGS)

Bioinformatics Support by Desiree Tillo, Bioinformatics Consultant


- Data delivery from NGS instruments

230217_VH01090_53_AACKM5HM5 Data Available

 mcintoshc@helix.nih.gov
To NCI LEC DNA Core

 Reply  Reply All  Forward  

Mon 2/20/2023 11:30 AM

 You forwarded this message on 2/20/2023 12:47 PM.
We removed extra line breaks from this message.

Content-Type: text/plain; charset="us-ascii"
MIME-Version: 1.0
Content-Transfer-Encoding: 7bit

Dear Core User,

Your NextSeq run on 2023-02-17 in data directory "230217_VH01090_53_AACKM5HM5" with experiment name "David Takeda Cell_Plasma ChIPSeq CS_RX 021723" completed successfully.

The data is accessible from the following URL:

[http://cleversafetest.nci.nih.gov/SEQ37V/NEXTSEQ2000/230217_VH01090_53_AACKM5HM5.tar?
Signature=kF2plqkRepHF8Xw%2F10S%2BT7v6%2BXg%3D&Expires=1678119203&AWSAccessKeyId=l6jcnKaAlMCbhfDQpECs](http://cleversafetest.nci.nih.gov/SEQ37V/NEXTSEQ2000/230217_VH01090_53_AACKM5HM5.tar?Signature=kF2plqkRepHF8Xw%2F10S%2BT7v6%2BXg%3D&Expires=1678119203&AWSAccessKeyId=l6jcnKaAlMCbhfDQpECs)

The QC file is accessible from the following URL:

[http://cleversafetest.nci.nih.gov/SEQ37V/NEXTSEQ2000/230217_VH01090_53_AACKM5HM5_qc.zip?
Signature=1XlzSyA7dClixg0jXOVH1UjT%2Ffc%3D&Expires=1678119203&AWSAccessKeyId=l6jcnKaAlMCbhfDQpECs](http://cleversafetest.nci.nih.gov/SEQ37V/NEXTSEQ2000/230217_VH01090_53_AACKM5HM5_qc.zip?Signature=1XlzSyA7dClixg0jXOVH1UjT%2Ffc%3D&Expires=1678119203&AWSAccessKeyId=l6jcnKaAlMCbhfDQpECs)

Notes on Obtaining your data ##### There are two ways to obtain your data files from an NIH Network attached computer:

- 1) Use the url link above.
- 2) Cut and paste the following command in a terminal window:

```
wget "http://cleversafetest.nci.nih.gov/SEQ37V/NEXTSEQ2000/230217_VH01090_53_AACKM5HM5.tar?  
Signature=kF2plqkRepHF8Xw%2F10S%2BT7v6%2BXg%3D&Expires=1678119203&AWSAccessKeyId=l6jcnKaAlMCbhfDQpECs" -O  
230217_VH01090_53_AACKM5HM5_Data.tar
```

Obtaining your QC data ##### There are two ways to

Service : Next Generation Sequencing (NGS)

Bioinformatics Support by Desiree Tillo, Bioinformatics Consultant

- Data delivery from NGS instruments
- Experimental design consultations
- Interpretation of initial sequencing QC
- Provides advice on analysis resources for a variety of genomics applications (ChIP-seq, RNA-seq, ATAC-seq, variant analysis, GeoMX DSP transcriptome atlas)
- Long-term research projects on a collaborative basis in certain cases
 - development of custom workflows/analyses
 - project and schedule dependent

Office of Science & Technology Resources (OSTR) Subsidy

- CCR investigators are eligible for NGS subsidy through OSTR from services from 3 NCI CCR Cores
 - Genomics Core (Bldg.41)
 - Single Cell Analysis Facility (SCAF) (Bldg.41)
 - Sequencing Facility (Fredericks)

Subsidies provided to PIs who exceed \$15,000 in cumulative sequencing costs during the fiscal year. After the initial \$15,000 threshold, cumulative sequencing costs will be subsidized at 50% up to a total of \$15,000 subsidy/year/PI

- CCR investigators are also eligible for Supplemental Technology Award Review System (STARS)
 - Pilot projects to get support for consumables

Thank You

We look forward to
**'Empowering your research
with services from the
CCR Genomics Core'**

ncilecdnacre@mail.nih.gov