

NCI Flow Cytometry: Other

e.g. fluid, aspirate (FNA), tissue

⌘ Before Collection Day:

- Request testing schedule:
 - Yuan, Constance yuanc@mail.nih.gov
 - Wang, Hao-Wei hao-wei.wang@nih.gov

CRIS order: **Flow Cytometry, Other**

⌘ On Collection Day: DELIVERY Cutoff 1:00 p.m.

Select specimen collection container by anticipated volume:

< 10 mL: Retrieve blue cap tube (3SWN Procedures unit)

10 - 300mL: Use plain collection container

- In SoftID, print specimen label

⌘ Transfer **< 10mL specimen** into CSF Media tube (blue cap tube);

⌘ Transfer **> 10mL specimen** into plain collection container

STAT Transport (dry) to Room 3S-240 / 3S-238



Where can I find FLOW information?



⌘ Flow Cytometry webpage:

<https://ccrod.cancer.gov/confluence/display/CCRLP/Flow+Cytometry>

Acting Supervisor: Truc Ho

Rm: 3S241 Office: 301-480-8834